



Instructions for Use

Membrane Oxygenator

OxyPrime ECMO

Adult

SUMMARY

1. DESCRIPTION	3
1.1. GRAPHIC IMAGES	4
1.2. TECHNICAL SPECIFICATIONS	5
2. PRESENTATION	6
3. INTENDED PURPOSE	6
4. INDICATION FOR USE	6
5. CONTRAINDICATIONS	6
6. CLINICAL BENEFITS	7
7. ADVERSE EFFECTS	7
8. INSTRUCTION FOR USE	7
8.1 SET-UP	7
8.2 CIRCUIT ASSEMBLY	7
8.3 FILLING WITH PRIME	8
8.4 START OF ECMO	10
8.5 ECMO CONTROL	10
8.6 END OF ECMO	10
8.7 REPLACING THE OXYGENATOR	10
9. WARNINGS	11
10. PRECAUTIONS	12
11. STORAGE AND TRANSPORT	14
12. PRODUCT DISPOSAL	14
12.1 SAFE DISPOSAL PROCEDURE	14
13. PERFORMANCE DATA	15
14. SHELF-LIFE	16
15. WARRANTY	17
16. MANUFACTURER INFORMATION	17
17. SYMBOLS DESCRIPTION	18
18. LABEL TEMPLATE	19

READ INSTRUCTIONS BEFORE USE

Model: Membrane Oxygenator Oxyprime ECMO ADULT (code 617320)

1. DESCRIPTION

The Membrane Oxygenator Oxyprime ECMO Adult consists of two coupled chambers (Gas Exchanger and Heat Exchanger), being indicated for long-term procedure. The internal space of the gas exchanger is filled with a Polymethylpentene (PMP) membrane, formed by capillaries (fibers) with a highly porous surface. The high porosity of the capillaries, associated with a uniform distribution of pores, provides the membrane with high gas transfer (diffusion) performance, desirable to meet the fundamental principle of an oxygenator: efficiency in transferring O₂ to the blood and in removing carbon dioxide (CO₂).

The internal space of the heat exchanger is filled with capillaries of thermoplastic polymer (polyester). This thermoplastic composition, combined with the characteristics of the production process, results in non-porous capillaries, totally impermeable to blood and water, and with high heat transfer capacity.

All the Oxygenator surfaces in contact with the blood are treated with the biocompatible coating, named BRCoating, composed of water produced by reverse osmosis, Polyethylene glycol 300 and recombinant human albumin*. This solution coats the inner surface of the product, reducing the interaction of the blood with it.

*Recombinant human albumin it is not derived from Animal or Human sources.

1.1. Graphic images




Figure 1 - Membrane Oxygenator Oxyprime ECMO Adult




Figure 2 - Dedicated Support (Holder) – Sold Separately

1.2. Technical Specifications

Table 1 - Technical Specifications

Material	
Housing/connector Material	Polycarbonate
Protective covers	Polyethylene
Coating	BRCoating
Performance	
Maximum Blood flow rate	7.0 L/min
Blood inlet maximum pressure	100 kPa (760 mmHg)
Priming volume (Static)	270 ml
Gas flow: blood flow ratio	0.5:1 - 1:1 - 2:1
Gas Exchanger	
Fiber Material	Polymethylpentene
Surface Area	2.0 m ²
Temperature connector	YSI
Blood outlet Direction: Right Side	3/8" (9.5 mm)
Gas inlet Connector	1/4" (6.4 mm)
Recirculation Connector	1/4" (6.4 mm)
Arterial Sample Connector	Luer lock
Removal air connector	Luer lock
Heat Exchanger	
Fiber material	Polyester
Surface Area	0.35 m ²
Water Maximum temperature	42°C
Water maximum pressure	32 PSI (1660 mmHg)
Blood inlet	3/8" (9.5 mm)
Water connector	1/2" (12.7 mm) Hansen quick connect fittings

Disclosable information: Information on blood cell damage, data tolerances and leaching.

If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

2. PRESENTATION

The Membrane Oxygenator Oxyprime ECMO Adult is sterile (sterilized by Ethylene Oxide - ETO) and non-pyrogenic. In its primary packaging, the product is packed in Tyvek paper envelope and polyethylene film. In its secondary packaging, the Membrane Oxygenator Oxyprime ECMO Adult is packed in a cardboard box. The instruction manual accompanies the product, in the secondary packaging.

3. INTENDED PURPOSE

The Membrane Oxygenator Oxyprime ECMO Adult is intended to be used in an extracorporeal perfusion circuit to oxygenate and remove carbon dioxide from the blood during extracorporeal membrane oxygenation (ECMO) technique in extracorporeal life support (ECLS) procedure. Braile does not specify the maximum period of use of the Membrane Oxygenator Oxyprime ECMO Adult, nor the ideal time to change it. The decision on whether to replace the product rests with the medical team, which should observe signs such as a decrease in the blood oxygenation rate and a significant increase in the pressure gradient.

4. INDICATION FOR USE

The Membrane Oxygenator Oxyprime ECMO Adult is indicated for the extracorporeal membrane oxygenation (ECMO) technique in the extracorporeal life support (ECLS) procedure in patients with cardiac and/or respiratory failure.

The Membrane Oxygenator Oxyprime ECMO Adult is for use with patients when required blood flow rate (cardiac output) will not exceed 7L/min.

This device shall only be used by properly trained and qualified personnel.

5. CONTRAINDICATIONS

The product is contraindicated to perform other functions that are not determined by the indication for use, or when, in the opinion of the medical team, its use is contrary to the best benefits for the patient.

6. CLINICAL BENEFITS

Use of this medical device may cause: Pressure drop, Hemolysis, Blood Loss, Inflammatory reaction, Embolism, Hypoxia, Infection Organ and tissue failure; Septic Shock, Kidney Failure, Heart attack, cell death, Permanent damage; Death.

7. ADVERSE EFFECTS

Use of this medical device may cause: Blood Loss, Hypoxia, Inflammatory reaction, organ and tissue failure; permanent damage; death.

8. INSTRUCTION FOR USE

It is recommended to use the Braile Biomédica Membrane Oxygenator Oxyprime ECMO Adult with devices from the Braile Biomédica brand, listed in Table 2.

Table 2 - Compatible devices (sold separately)

Description	Code
Tubing Set ECMO Adult	617327
Centrifugal Pump Safira Centriflux with BRCoating	616322

* Braile Biomédica cannot guarantee the compatibility between products from other manufacturers and, consequently, their safety and effectiveness.

8.1 Set-Up

- a. **Support (Holder):** Turn the knob counterclockwise to free the support Fixing clamp and be able to position it on the rod freely. Set the holder at the appropriate height to prevent the lines from kinking, and turn the handle clockwise to fix the support. The oxygenator support must always be fixed below patient's level;
- b. **Oxygenator:** Open the sterile packaging using aseptic techniques.
- c. Remove the oxygenator and attach it to its support, by fit the hole in the oxygenator base into the support base fitting. Subsequently, fit the upper part of the oxygenator (gas exchanger) into the movable claws of the support, make sure that the oxygenator is well fixed and secure in the support.

8.2 Circuit Assembly

- a. Connect the ECMO arterial line from the tube set to the arterial outlet of the oxygenator.
- b. Connect the ECMO venous line from the tube set to the blood inlet of the centrifugal pump.

- c. Connect the centrifugal line of the tube set between the blood outlet of the centrifugal pump and the blood inlet of the oxygenator.
- d. Remove the centrifuge from the device coupling and hang it up using the tubes.
- e. Connect the oxygen line of the tube set between the gas inlet of the oxygenator and the gas mixer of the ECMO device.
- f. Connect the 3-way stopcocks of the tube set to the luer connectors of the prime connectors 1 and 2 of the venous line of the tube set.
- g. Connect the serum lines in the prime connectors 1 and 2, through the 3-way stopcock, and then connect to a serum bag with two openings for infusion sets.
- h. Connect the oxygenator recirculation line between the bubble-trap connector and recirculation.
- i. Connect an MM extender of the oxygenator between the stopcock of the bubble-trap connector and the stopcock of the prime connector 2.
- j. Place a tube clamp between the prime connectors 1 and 2. Open the flow regulating clamps of the serum lines 1 and 2.
- k. According to the team routine for pressure measurements, place 3-way stopcocks in the circuit luer connectors.
- l. Connect the water lines from the ECMO device to water connectors in the oxygenator.

Note 1: Read the instructions on the manual for the Centrifugal Pump Safira Centriflux with BRCoating and Tubing Set ECMO Adult for complete instruction on circuit assembly.

8.3 Filling with Prime

- a. Allow serum to flow through serum line 1, filling the centrifugal pump and oxygenator by gravity. Check that the oxygenator recirculation and bubble-trap lines are open.
- b. With the centrifugal pump and oxygenator filled, close the oxygenator recirculation and bubble-trap lines.
- c. Place the centrifugal pump in the ECMO device coupling and start rotation.
- d. Once the entire circuit is filled and free of bubbles, close serum lines 1 and 2 and remove the clamp between them.
- e. Stop rotation of the console; disconnect the serum lines from the prime connectors 1 and 2. Close the luer with a sterile luer cover. Make sure there are no air bubbles in the body of the luer.

8.4 Start of ECMO

- a. Remove the clamps from the venous line.
- b. Increase the ECMO device rotation speed to 1700rpm.
- c. Remove the clamp from the arterial line simultaneously with the increase in rotation.
- d. Adjust the blood flow according to the patient's need.
- e. Start the gas flow and adjust to the desired gas flow to blood flow ratio (0.5: 1 - 1: 1 - 2: 1).

8.5 ECMO control

- Monitor the performance of the oxygenator by periodic blood gas analysis.
- Monitor the patient clotting status by periodically analyzing TCA.
- Monitor the pressure gradient of the oxygenator periodically.
- Control the blood temperature by regulating water temperature, ECMO device heater, which circulates in the oxygenator heat exchanger.

Problems and solutions	
low pO ₂	Increase O ₂ %
high pO ₂	Reduce O ₂ %
low pCO ₂	Increase gas flow
high pCO ₂	Reduce gas flow

8.6 End of ECMO

- a. Decrease the rotation of the ECMO device, and consequently blood flow.
- b. Clamp the arterial and venous lines simultaneously.
- c. Close the gas flow.
- d. Stop the rotation of the ECMO device.
- e. Turn off heating / water flow of the ECMO device.

8.7 Replacing the Oxygenator

- a. Keep a spare oxygenator.
- b. Clamp the arterial and venous lines simultaneously.
- c. Stop rotation and water flow in the ECMO device.
- d. Close gas flow.
- e. Place tube clamps close to the oxygenator blood outlet and blood inlet connectors, close the recirculation line.

- f. Remove the oxygenator to be replaced from the holder and replace with the new one.
- g. Remove the arterial line from the oxygenator to be replaced and connect it to the new one.
- h. Remove the venous line from the oxygenator to be replaced and connect to the new one.
- i. Place the recirculation line in the new oxygenator and keep it open.
- j. Open the clamp next to the blood inlet connector and allow the solution to flow until it is completely filled, through the prime inlet connector.
- k. Make sure there are no bubbles in the circuit.
- l. Start ECMO according to the ECMO start instructions (**Item 8.4**).

9. WARNINGS

- Read the instructions before use;
- Single use product. Therefore, do not re-sterilize or reuse the product. There are no data to support the sterility, pyrogenicity and functionality of devices after reprocessing;
- Reprocessing prohibited. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device (Leaks, low performance and membrane rupture);
- Do not use non-sterile products - they must be sterile before use;
- Do not use the product after the expiry date;
- Do not use the product if it is damaged or the packaging is broken;
- Before use, the product must be checked for proper operation and to ensure the model is suitable for its intended purpose;
- Product intended exclusively for medical use;
- The product should only be used by people specially trained on ECMO techniques.
- The procedure and the correct application of the product are the responsibility of the doctor who performs the treatment;
- The product must not be used under conditions that contradict its technical specifications, under penalty of loss of performance or accidents. It is up to the person responsible for using the product to check in advance whether the technical specifications meet the needs of the procedure to be performed;
- Always use aseptic techniques when handling the product after opening;
- Avoid any impacts during transport or use;
- Do not use the product in case of doubt;
- The product should be used immediately after opening the package. After use,

dispose of safely as medical waste according to health institution policies. The product is biohazardous as it is contaminated with blood;

- Pressure at the blood inlet of the oxygenator should not exceed 100kPa (760 mmHg). Pressure greater than 100kPa may cause leaks or damage to the device;
- Water pressure at the heat exchanger inlet should not exceed 220kPa (1660mmHg). Pressure greater than 220kPa may cause leaks or damage to the device;
- When using the centrifugal pump on the arterial line, clamp the arterial line distal to the oxygenator (the patient's side) before stopping the pump. Improper clamping may cause back-flow or migration of gaseous emboli into the blood side;

10. PRECAUTIONS

- It is recommended that you check the product carefully before use. The operation of the product packaging may be affected by transport. It is not possible to guarantee full operation of the product when transport damage occurs. If the product has been dropped or crushed, it cannot be used and must be replaced with a new product;
- Always keep a spare unit of the product during the procedure;
- The manufacturer is not responsible for damage resulting from lack of experience, misuse or failure to follow the instructions for use;
- During the procedure, blood anticoagulation protocol should be adopted;
- Protamine should not be used in patients before or during extracorporeal circulation;
- It is necessary to provide correct heparinization and accurate monitoring of the anticoagulation status before, during or after extracorporeal circulation;
- Check all circuit connections in advance to prevent any leaks;
- Remove all air from the extracorporeal circuit before the start of ECMO using the circuit prime procedure;
- It is recommended to start the procedure immediately after filling with prime in order to minimize the period of microbiological contamination, precipitation of the filling solution;
- Do not make changes to the product;
- Do not use the product if not all connectors are protected with a protective cover. The protective covers must be retained immediately before use;
- Regularly monitor the patient to identify side effects of cardiopulmonary bypass (e.g., infections, hemolysis, post-perfusion syndrome and organ damage);
- Regularly monitor the patient clotting status (e.g., activated clotting time and partial thromboplastin time). The anticoagulation control protocol is the responsibility of the responsible user;

- Broken and / or leaking connectors must be replaced;
- Tube metal clamps should be available;
- During transport of patient, avoid leaving the tubes hanging, excessive forces, bumps and kinks;
- To collect samples, avoid creating bubbles in the circuit;
- Monitor collection and recirculation lines to prevent clots;
- To ensure optimal operation of the oxygenator, additional devices or connections are required. Read the instructions for use carefully before using the device, ensuring that all necessary parts are available before use;
- Avoid bumping into the oxygenator during the prime procedure;
- During infusion, the blood pressure in the oxygenator gas exchanger must always be kept higher than the gas pressure to prevent formation of micro-bubbles in blood (gas embolism);
- The gas outlets of the oxygenator must never be blocked, as this causes an increase in gas pressure in the gas exchanger, which can cause gas embolism;
- The maximum blood flow of the oxygenator model used must not be exceeded;
- Blood flow should never be reduced suddenly. The sudden reduction or sudden stop of blood flow causes a decrease in blood pressure in the gas exchanger, which can cause gas embolism. The reduction or stop of blood flow must always be done gradually and accompanied by reduction in gas flow, with continuous monitoring of gas and blood pressures in the gas exchanger;
- Whenever the Arterial Pump is stopped (zero flow), the gas flow must be closed (zero flow), under the risk of gas embolism;
- Monitor oxygenation and carbon dioxide removal throughout the procedure through periodic gas analysis;
- Monitor the suction pressure in the venous line, check if luer connections and 3-way valves are closed to avoid embolism in the line;
- Monitor oxygenator pressure gradient (pressure before and after the oxygenator);
- Monitor gas exchange performance in gas analysis;
- Always position the oxygenator below the patient level;
- Never infuse drugs directly into the oxygenator connectors;
- Do not use solvents such as alcohol, ether, acetone, liquid anesthetics by inhalation (for example: Halothane, Enflurane, etc.), as these can damage the product;
- In case of need to replace the Oxygenator during procedure, proceed according to related instructions contained in the Product Replacement Instructions.

11. STORAGE AND TRANSPORT

The product in its original packaging must be stored in a clean, dry, weather-free environment and at temperatures above -10°C and below 40°C. For its transportation, it should never be refrigerated, frozen or exposed to temperatures above 40°C.

12. PRODUCT DISPOSAL

After use, follow the **safe disposal procedure, item 12.1**, to safely dispose of the device. Product and packaging disposal must comply with the hospital and / or local government policy. The product must be disposed of only in hospital / infectious waste. The product is biohazardous as it is contaminated with blood.

12.1 Safe Disposal Procedure

Warning: Wear gloves to carry out this process, there is a risk of contamination.



Note 2: It is not necessary to disconnect the bloodlines tubes.



- a. Detach the centrifugal pump and set aside on a secure surface.
- b. Position a surgical field below the oxygenator to avoid wetting the floor, clamp the water lines near the oxygenator's water inlet/outlet connectors, and, disconnect the lines, set the pipes aside on a safe surface.

Note 3: For the next steps, use the disposal waste for contaminated products.

- c. Remove the oxygenator from the support (holder) and discard it carefully.
- d. Discard the centrifugal pump and pump line.
- e. Discard the remaining lines.

13. PERFORMANCE DATA





14. SHELF-LIFE

The Membrane Oxygenator Oxyprime ECMO ADULT is valid for 2 years after the sterilization date, considered sterile only if the package has not been violated. Check the validity date on the label of the product package.

15. WARRANTY

Braile Biomédica Indústria, Comércio e Representações Ltda, certifies that the product has been manufactured with due care, with good quality materials, and is warranted for manufacturing defects pursuant to the Brazilian Consumer Protection Code. It further clarifies that if the product is violated by unauthorized or unqualified personnel, the warranty will be void. The manufacturer states that any other warranty expressly or implicitly granted by others shall not be valid. The use of the product is restricted to doctors or qualified and trained technicians under the supervision of the responsible physician, who are familiar with the procedures for using the product, being their responsibility said use. The manufacturer will not be responsible for any damage resulting from product misuse.

Braile Biomédica Ltda also states that the right to complain about the product will be governed by the provisions contained in Article 26 of the Brazilian Consumer Protection Code (Law 8.078 / 90). No agent, employee, representative, or distributor of the Braile Biomédica Ltda has the authority to change or amend what is described herein, assume or bind the Braile Biomédica Ltda to any other liability or warranty in connection with this product.

This product was developed and manufactured based on legal and regulatory requirements, and under strict quality control, aiming at managing potential risks of its use, in order to minimize them within safe limits. However, Braile Biomédica cannot ensure that the product is totally free from the likelihood of failure, so that its use must respect all precautions, restrictions, warnings, special care and relative instructions for use, as well as must be constantly and carefully monitored.

16. MANUFACTOR INFORMATION

Manufactured by:

Braile Biomédica Indústria Comércio e Representações Ltda

CNPJ N°: 52.828.936/0001-09

Av. Presidente Juscelino Kubitschek de Oliveira, 1505 – Jardim Tarraf I

CEP: 15091-450

São José do Rio Preto – São Paulo

Imported and distributed by:

Mars Medical

























Landhausstrasse 46

70190, Stuttgart – Germany

Tel: +49 1751938653

E-Mail: info@marsmedical.de

17. SYMBOLS DESCRIPTION

	Manufacturer		Use-by date
	Catalog Number		Sterile
	Serial Number		Sterilized using ethylene oxide
	Batch Code		Do not use if packaging is damaged
	Date/ Country Manufacture		Fragile, handle with care
	Non-pyrogenic		This way up
	Consult Instructions for Use or electronic instructions for use		Stacking limit
	Temperature limit		Recyclable
	Keep dry		Medical Device
	Keep away from sunlight		Do not re-use
	Unique Device Identifier		Do not re-sterilize
	CE Mark		Authorized representative in the European Union

18. LABEL TEMPLATE

	<p>Braile Biomédica Indústria Comércio e Representações Ltda. Av. Presidente Juscelino Kubitschek de Oliveira, 1505 - Jardim Tarraf I CEP 15091-450 - São José do Rio Preto - SP - Brasil Fone: +55 17 2136-7000 SAC: 0800 707 20 50 CNPJ: 52.828.936/0001-09 www.braile.com.br</p>							
<table border="1" style="border-collapse: collapse;"> <tr> <td style="padding: 2px;">EC</td> <td style="padding: 2px;">REP</td> </tr> </table>	EC	REP	<p>Mars Medical Landhausstrasse 46 70190, Stuttgart – Germany Tel: +49 1751938653 E-Mail: info@marsmedical.de</p>					
EC	REP							
<p>MEMBRANE OXYGENATOR OXYPRIME ECMO ADULT</p>								
<p>Description (Model): 617320 Package Content: 01 unit of Membrane Oxygenator Oxyprime ECMO Adult / 01 Instruction for Use Contains 1 unit destined to medical application</p>								
<p>ANVISA Registration nº: XXXXX Technical Name: Oxygenator Technician in Charge: Vladimir D. A. Ramirez – CRF-SP 09010</p>								
<p>Read instructions before use Single use product / Prohibited reprocessing</p>								
<p>Indications, Precautions, Warnings, Instructions of Use, Preservation Cares*. * Please refer to leaflet insert "INSTRUCTIONS OF USE".</p>								
<table border="1" style="border-collapse: collapse;"> <tr> <td style="padding: 2px;">REF</td> <td style="padding: 2px;">XXXXXX</td> </tr> </table>	REF	XXXXXX	<table border="1" style="border-collapse: collapse;"> <tr> <td style="padding: 2px;">SN</td> <td style="padding: 2px;">XXXXXX</td> </tr> </table>	SN	XXXXXX	<table border="1" style="border-collapse: collapse;"> <tr> <td style="padding: 2px;">LOT</td> <td style="padding: 2px;">XXXXXX</td> </tr> </table>	LOT	XXXXXX
REF	XXXXXX							
SN	XXXXXX							
LOT	XXXXXX							
<table border="1" style="border-collapse: collapse;"> <tr> <td style="padding: 2px;">  </td> <td style="padding: 2px;">YYYY/MM/DD</td> </tr> </table>		YYYY/MM/DD	<table border="1" style="border-collapse: collapse;"> <tr> <td style="padding: 2px;">UDI</td> </tr> </table>	UDI				
	YYYY/MM/DD							
UDI								
<table border="1" style="border-collapse: collapse;"> <tr> <td style="padding: 2px;">  </td> <td style="padding: 2px;">YYYY/MM/DD</td> </tr> </table>		YYYY/MM/DD						
	YYYY/MM/DD							
<table border="1" style="border-collapse: collapse;"> <tr> <td style="padding: 2px;">STERILE</td> </tr> </table>	STERILE	<table border="1" style="border-collapse: collapse;"> <tr> <td style="padding: 2px;">STERILE</td> <td style="padding: 2px;">EO</td> </tr> </table>	STERILE	EO				
STERILE								
STERILE	EO							
			<table border="1" style="border-collapse: collapse;"> <tr> <td style="padding: 2px;">MD</td> </tr> </table>	MD	<p>braile.com.br/en/ produtos/cemark products</p>			
MD								

The external package has the following informations

The product must be stored and transported in its original packaging, in a clean, dry, weather-free environment and at a controlled temperature.

