

WIRED VENOUS CANNULA



Thin Wall
Transparent and Flexible Body
Easy and Accurate Insertion
Kink Resistant

Description: The Wired Venous Cannula is made of overlapped layers of transparent and flexible plastic polyvinyl chloride-based resin (PVC). The tip may have a straight or angled shape, with one or two steps of drainage; the body has a cylindrical or oval shape and it is reinforced by a stainless steel wire, helical shape (anti-kink), positioned among the plastic layers, without any blood contact. The metallic tip has a beveled shape with a side hole. Each package contains a sterile and pyrogen-free unit.

Indication: The Wired Venous Cannula is intended for venous drainage in surgical procedures performed with cardiopulmonary bypass.

Description	Diameter	Length	Connection	Connector	Code	Situation
Wired Venous Cannula - Single -Stage Straight Tip	12 FR	31 cm	1/4"	—	162274	Standard Product
Wired Venous Cannula - Single -Stage Straight Tip	14FR	31 cm	1/4"	—	162283	Standard Product
Wired Venous Cannula - Single -Stage Straight Tip	16 FR	31 cm	1/4"	—	162292	Standard Product
Wired Venous Cannula - Single -Stage Straight Tip	18 FR	36 cm	3/8"	—	162308	Standard Product
Wired Venous Cannula - Single -Stage Straight Tip	20 FR	36 cm	3/8"	—	162317	Standard Product
Wired Venous Cannula - Single -Stage Straight Tip	22 FR	36 cm	3/8"	—	162326	Standard Product
Wired Venous Cannula - Single -Stage Straight Tip	24 FR	36 cm	3/8"	—	162335	Standard Product
Wired Venous Cannula - Single -Stage Straight Tip	28 FR	40 cm	3/8"	—	157076	Standard Product
Wired Venous Cannula - Single -Stage Straight Tip	30 FR	40 cm	3/8"	—	139959	Standard Product
Wired Venous Cannula - Single -Stage Straight Tip	32 FR	40 cm	3/8"	—	139968	Standard Product
Wired Venous Cannula - Single -Stage Straight Tip	34 FR	40 cm	3/8"	—	139977	Standard Product
Wired Venous Cannula - Single -Stage Straight Tip	36 FR	40 cm	1/2"	—	162344	Standard Product
Wired Venous Cannula - Single -Stage Tip 90°	12 FR	27 cm	1/4"	—	224536	Made to Order
Wired Venous Cannula - Single -Stage Tip 90°	14 FR	27 cm	1/4"	—	224545	Made to Order
Wired Venous Cannula - Single -Stage Tip 90°	16 FR	26 cm	1/4"	—	224554	Made to Order
Wired Venous Cannula - Single -Stage Tip 90°	18 FR	31 cm	3/8"	—	224563	Made to Order
Wired Venous Cannula - Single -Stage Tip 90°	20 FR	31 cm	3/8"	—	224572	Made to Order
Wired Venous Cannula - Single -Stage Tip 90°	22 FR	31 cm	3/8"	—	224581	Made to Order
Wired Venous Cannula - Single -Stage Tip 90°	24 FR	31 cm	3/8"	—	224606	Made to Order
Wired Venous Cannula - Single -Stage Tip 90°	28 FR	34 cm	3/8"	—	224624	Standard Product
Wired Venous Cannula - Single -Stage Tip 90°	30 FR	34 cm	3/8"	—	224633	Made to Order
Wired Venous Cannula - Single -Stage Tip 90°	32 FR	34 cm	3/8"	—	224642	Standard Product
Wired Venous Cannula - Single -Stage Tip 90°	34 FR	35 cm	3/8"	—	224660	Made to Order
Wired Venous Cannula - Single -Stage Tip 90°	36 FR	34 cm	1/2"	—	224679	Made to Order
Wired Venous Cannula - Double - Stage	28/36 FR	39 cm	3/8"	—	235477	Made to Order
Wired Venous Cannula - Double - Stage	28/36 FR	39 cm	1/2"	—	242617	Made to Order
Wired Venous Cannula - Double - Stage	29/37 FR	38 cm	1/2"	—	242547	Made to Order
Wired Venous Cannula - Double - Stage	29/29 FR	39 cm	3/8"	—	242556	Made to Order
Wired Venous Cannula - Double -Stage	32/40 FR	31 cm	1/2"	—	235459	Standard Product
Wired Venous Cannula - Double -Stage	34/46 FR	31 cm	1/2"	—	235440	Standard Product
Wired Venous Cannula - Double -Stage	36/46 FR	31 cm	1/2"	—	235468	Standard Product
Wired Venous Cannula - Double -Stage	36/51 FR	39 cm	1/2"	—	270830	Standard Product
Wired Venous Cannula - Double -Stage	36/50 FR	45 cm	1/2"	1/2" x 1/2"	131441	Made to Order
Wired Venous Cannula - Double -Stage Oval	32/40 FR	31 cm	1/2"	—	242565	Made to Order
Wired Venous Cannula - Double -Stage Oval	36/46 FR	32 cm	1/2"	—	235486	Standard Product
Wired Venous Cannula - Double -Stage Oval	36/51 FR	39 cm	1/2"	—	270821	Made to Order
Wired Venous Cannula - Double -Stage Oval	36/50 FR	45 cm	1/2"	1/2" x 1/2"	139986	Made to Order

*±1cm

Código Catálogo 346643
Rev. 02

Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
Zipe Code 15091-450 • São José do Rio Preto • SP • Brazil
Phone #55 17 2136 7000 • Fax #55 17 2136 7040
Federal Tax ID: 52.828.936/0001-09
Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010
ANVISA Reg. nº 10159030092
www.braille.com.br

