

Cannula for Antegrade Cardioplegia



www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA

Description: The Cannula for Antegrade Cardioplegia comprises a PVC tube with a vent line and flow regulator, metallic introducer needle and tip with side drilling and flange (suture ring) for aiding the attachment. A female luer lock adapter in its end permits the connection with the cardioplegia system, and a 1/4" straight connector permits the coupling of ventline. Each package contains a sterile and pyrogen-free unit.

Indication: The Cannula for Antegrade Cardioplegia is intended for infusion of cardioplegic solution, using antegrade route in procedures with cardiopulmonary bypass.

Code

182643 - Cannula for Antegrade Cardioplegia

Technical Specifications

Cannula Material	PVC
Diameter of Cannula Body	12 Fr
Cannula Length	20 cm
Tip Diameter	9 Fr
Needle Material	Aço Inox
Needle Specifications	14 G
Cannula Adapter	<i>luer lock</i> fêmea
Vent Line Connector	1/4"
Vent Line Length	10,5 cm

Catalogue Code 600250
Rev



Braille Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
CEP 15091-450 • São José do Rio Preto • SP • Brasil
Telefone 17 2136 7000 • Fax 17 2136 7040
Other countries • Phone # (55) 17 2136 7000 • Fax 17 2136 7040
CNPJ: 52.828.936/0001-09
Responsável Técnico: Vladimir D. A. Ramirez • CRF/SP 9010
Registro ANVISA nº 10159030082
www.braille.com.br

