


Balloon Catheter



www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA



Compliant balloon catheter with 2Fr profile and inflation diameter up to 45mm.



| Volume (ml) | Diameter (mm) |
|-------------|---------------|
| 10 | 20 |
| 18 | 25 |
| 30 | 35 |
| 40 | 40 |
| 60 | 45 |

* Refer to availability of other dimensions.

Intended for accommodating stent-graft and temporary vessel occlusion.



Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
 Zipe Code 15091-450 • São José do Rio Preto • SP • Brazil
 Phone #55 17 2136 7000 • Fax #55 17 2136 7040
 Federal Tax ID: 52.828.936/0001-09
 Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010
 ANVISA Reg. nº 10159030066
www.braile.com.br

BRAILE[®]
 BIOMÉDICA