

Mediastinal Drainage System



www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA

Description: The Mediastinal Drainage System is made of a graduated reservoir with a lid, draining line, vacuum line and catheter. Made in rigid transparent PVC with different volume capacities: 2000, 1000, 500 or 200 ml, still contains internal mariott tube, flow regulator and fastening cable.

The air and a body fluids are removed from the patient through the multiperforated catheter and drainage line, being deposited inside the bottle.

Indication: The Mediastinal Drainage System is indicated to remove air and liquid which can accumulate after surgical acts or certain clinical situations os hemothorax, pneumothorax and empyema.

Technical Specifications

System reservoir for Mediastinal Drainage

Reservoir raw material	PVC
Drainage line raw material	PVC
Vaccum line raw material	PVC
Internal Mariott tube raw material	Polyethylene
Connectors raw material	Polycarbonate
Reservoir volume	200, 500, 1000 e 2000

System catheter for Mediastinal Drainage

Catheter diameter	20, 22, 26 32, 36, 38 Fr
Raw material	PVC

Code

- 48237
Mediastinal Drainage System - 2000/32
- 38797
Mediastinal Drainage System - Adult 2000/36
- 62550
Mediastinal Drainage System - Infant 2000/26
- 101994
Mediastinal Drainage System - Infant 1000/26
- 21223
Mediastinal Drainage System - Adult 1000/38
- 34454
Mediastinal Drainage System - Adult 1000/38 SPB
- 117858
Mediastinal Drainage System - Adult 2000/38 PA
- 40527
Mediastinal Drainage System - Adult 2000/36 CIR
- 62567
Mediastinal Drainage System - Adult 2000/38
- 110529
Mediastinal Drainage System - Infant 200/2026
- 108384
Mediastinal Drainage Systeml - Infant 500/20
- 38340
Pressure Reducer 2000 ML



Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
Zipe Code 15091-450 • São José do Rio Preto • SP • Brazil
Phone #55 17 2136 7000 • Fax #55 17 2136 7040
Federal Tax ID: 52.828.936/0001-09
Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010
ANVISA Reg. nº 10159030083
www.braile.com.br

BRAILE
BIOMÉDICA