

Kit for Coronary Surgery without CEC Midcab - Kit



www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA

Description: The Kit for Coronary Surgery without Cardiopulmonary Bypass Midcab-Kit comprises a case which contains several disposable items used in the different techniques of this surgical procedure.

Indication: The Kit for Coronary Surgery without Cardiopulmonary Bypass Midcab-Kit is indicated for stabilization and occlusion of the coronary artery when performing anastomosis.

Code

103185
Kit for Coronary Surgery without CEC
Midcab - Kit

Technical Specifications

Components	Quantity (unit)
Cord no.	02
Mosquito Protector 20 cm Reperfusion	01
Instrument Flat Lace 10 mm Popcorn	02
Extensor MM	01
Silicone Tube	04
Thread Guide	02
Vessel Cannula	03
Bulldog	02
Jelco® nº 14	01
Blower CO2	01
Micro-Retractor of Coronary Artery	02
Silicone Thread 2.5/5 mm	02
Silicone Thread 1.25/2.5 mm	01
Shunt 1.5/18 mm	01
Shunt 1.75/19 mm	01
Shunt 2.0/20 mm	01
Occluder 150 mm	01
Occluder 175 mm	01
Occluder 200 mm	01



Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
Zipe Code 15091-450 • São José do Rio Preto • SP • Brazil
Phone #55 17 2136 7000 • Fax #55 17 2136 7040
Federal Tax ID: 52.828.936/0001-09
Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010
ANVISA Reg. nº 10159030038
www.braile.com.br

BRAILE[®]
BIOMÉDICA