

OxyPrime[®]


Membrane Oxygenator



- ▶ Low priming volume
- ▶ New *design*
- ▶ More compact
- ▶ High performance





www.braile.com.br

 **BRAILE[®]**
BIOMÉDICA



Techniques Especificacion

Material	Polycarbonate
Maximum blood flow	7,0 l/min
Total volume of priming	250 ml
Gas exchanger	
Fiber material	Polypropylene
Effective membrane area	1,95 m ²
Temperature sensor	YSI
Blood outlet connector	3/8"
Gas inlet connector	1/4"
Recirculation connector	luer lock
Arterial sample connector	luer lock
Air removal connector	luer lock
Heat exchanger	
Fiber material	Polyester
Effective membrane area	0,4 m ²
Maximum water temperature	42°C
Maximum water pressure	40 PSI
Blood inlet connector	3/8"
Water inlet/outlet connector	Hansen (1/2")



Commercial information

OxyPrime 7000

Code 610739

For orders or more information about any of our products, accessories or therapeutic solutions, contact your local sales representative.

Braile Biomédica Indústria, Comércio e Representações Ltda.
 Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
 Zipe Code 15091-450 • São José do Rio Preto • SP • Brazil
 Phone #55 17 2136 7000 • Fax #55 17 2136 7040
 Federal Tax ID: 52.828.936/0001-09
 Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010
 ANVISA Reg. n° 10159030098
 www.braille.com.br

BRAILE[®]
 BIOMÉDICA

