

Digital Triple Chronometer C3000

Simplicity
Accuracy
Reliability



The Digital Triple Chronometer c3000 Braile Biomedica comprises three high accuracy chronometer coupled to a single box and having independent operation. It works with a 9.6 V built-in battery which can be recharged and/or connected to the power grid, providing a continuous operation even when the local network power source is not available; microprocessed circuit for better reliability and clamp for attachment with two positions (vertical and horizontal). The Chronometer is intended for recording the time of the many surgical steps, such as: time of surgery, time of campling, time of cardioplegia, etc.

Code

353856

Digital Triple Chronometer



Technical Specifications

Accuracy	99,0% (25°C)
Display capacity	99:59:59 (h:m:s)
Duration of Battery	2 hours power grid or battery

Catalogue Code 205124
Rev. 01 | fev. 2015



Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
CEP 15091-450 • São José do Rio Preto • SP • Brasil
Telefone 17 2136 7000 • Fax 17 2136 7040
Other countries • Phone # (55) 17 2136 7000 • Fax 17 2136 7040
CNPJ: 52.828.936/0001-09
Responsável Técnico: Vladimir D. A. Ramirez • CRF/SP 9010

www.braille.com.br

 **BRAILE**[®]
BIOMÉDICA