

GUIDE WIRE VELOZ PLUS

HYDROPHILIC

Nitinol coated in hydrophilic polyurethane. Excellent radiopaque visualization.
Stiffness level: standard. Facilitated navigation. Safe positioning.
0.035 diameter (0.89 mm).



| DESCRIPTION | GUIDE WIRE VELOZ PLUS 0,035" x 260 cm (adaptable type) |
|--------------------|---|
| CODE | 607366 |
| COATING | Hydrophilic |
| MATERIAL OF CORE | Nitinol |
| FORMAT OF TIP | Adaptable |
| STIFFNESS OF TIP | Standard |
| DIAMETER OF WIRE | 0.035" |
| HARDNESS OF BODY | Standard |
| LENGTH OF THE WIRE | 260 cm |

ANVISA: nº 10159030100

GUIDE WIRE VELOZ

TEFLONIZED

Stainless steel coated in PTFE.
Extremely flexible and radiopaque distal tip end in J.
Stiffness levels: Standard or Extra Stiff.
Diameter of 0.035" (0.89 mm).

EXTRA STIFF

STANDARD

| DESCRIPTION | GUIDE WIRE VELOZ 0,035* x 260 cm (J3) Extra Stiff Teflon | GUIDE WIRE VELOZ 0,035* x 260 cm (J3) Standard Teflon |
|----------------|---|--|
| CODE | 606495 | 606500 |
| MODEL | Lunderquist | Seldinger |
| BODY RIGIDITY | Extra rigid | Standard |
| FORMAT OF TIP | J3 | J3 |
| COATING | Teflonized | Teflonized |
| FERRULE LENGTH | 5 cm | 4,5 cm |
| WIRE DIAMETER | 0.035" | 0.035" |
| WIRE LENGTH | 260 cm | 260 cm |

ANVISA: nº 10159030099

Braile Biomédica Indústria, Comércio e Representações Ltda.
Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
CEP 15091-450 • São José do Rio Preto • SP • Brasil
Telefone 17 2136 7000 • Fax 17 2136 7040
Other countries • Phone # (55) 17 2136 7000 • Fax 17 2136 7040
CNPJ: 52.828.936/0001-09
Responsável Técnico: Vladimir D. A. Ramirez • CRF/SP 9010



BRAILE
BIOMÉDICA