

Temperature Sensor



www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA

The temperature sensors are internally shielded thermoplastic cables, with a five tract connector with gold terminals and a small-sized PT-100 type pad, enabling faster measurements.

Technical Specifications

Scale	0 - 70°C
Connection	Redel 5 Vias
Accuracy	± 0,4°C
Length	3m
Sterilizable in Ethylene Oxide	

Code

117168 Temperature Sensor



Catalog Code 600329
Rev.



Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
Zip Code 15091-450 • São José do Rio Preto • SP • Brazil
Phone #55 17 2136 7000 • Fax #55 17 2136 7040
Federal Tax ID: 52.828.936/0001-09
Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010
ANVISA Reg. nº 10159030075
www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA