

Bovine Pericardium Patch Monocuspid




Experience
Quality
Performance

Description:

The Bovine Pericardium Patch Monocuspid is made with bovine pericardium with bovine pericardium cuspid, treated with glutaraldehyde and preserved in formaldehyde 4%. This treatment gives it proper characteristics of resistance, flexibility and lack of antigenicity.

Indication:

The Bovine Pericardium Patch Monocuspid is used in reconstructing the outflow tract of the right ventricle and pulmonary.

**TECHNICAL SPECIFICATIONS (cm)**

Size	A	B	ØC	Code
7	110	60	7	149000
9	110	60	9	144120
11	110	60	11	144139
13	110	60	13	144148
15	110	60	15	144157
17	110	60	17	144166
19	110	60	19	30049
21	110	60	21	30056

Code Catalog 600228
Rev. 01 - 07/2018

