

Organic Tubular Arterial Graft Valvulated Tube 200 Corrugated



Experience
Quality
Performance


 **BRAILE**[®]
BIOMÉDICA

Description:

The Organic Tubular Arterial Graft Valvulated Tube 200 Corrugated is made with a bovine pericardium molded in a special support, treated with glutaraldehyde and conserved in formaldehyde 4%. This treatment gives it proper characteristics of resistance, flexibility and lack of antigenicity.

Indication:

The Organic Tubular Arterial Graft Valvulated Tube 200 Corrugated is indicated for corrections of several cardiopathies, with aneurysms of the ascending aorta with or without valvular insufficiency, atrioventricular congenital hipoplasia, transposition of great arteries, intraventricular communication, annuloaortic ectasia, obstruction of the left ventricular outflow, stenosis and pulmonary atresia, among others.



TECHNICAL SPECIFICATIONS						
Size	A ± 10	ØB ± 0,5	ØC ± 1,0	D ± 0,5	Valve	Code
19	120	19	17	1,7	Ø19	11611
21	120	21	19	1,7	Ø21	11628
23	120	23	21	1,7	Ø23	11635

The size of the tube is specified by the valve Ø.

Code Catalog 600218
Rev.01 - 07/2018

