

# Organic Tubular Arterial Graft Valvulated Tube 200 Smooth



Experience  
Quality  
Performance


 **BRAILE**<sup>®</sup>  
BIOMÉDICA

**Description:**

The Organic Tubular Arterial Graft Valvulated Tube 200 Smooth is made of bovine pericardium treated with glutaraldehyde and conserved in formaldehyde 4%. This treatment gives it proper characteristics of resistance, flexibility and lack of antigenicity.

**Indication:**

The Organic Tubular Arterial Graft Valvulated Tube 200 Smooth is indicated for corrections of several cardiopathies, such as aneurysms of the ascending aorta with or without valvular insufficiency, atrioventricular congenital hipoplasia, transposition of major arteries, intraventricular communication, annuloaortic ectasia, obstruction of the left ventricular outflow tract, stenosis and pulmonary atresia, among others.



TECHNICAL SPECIFICATIONS (mm)						
Size	A(mm) ± 10	ØB(mm) ± 0,5	ØC(mm) ± 0,5	D(mm) ± 0,5	Valve	Code
19	170	19	18	1,7	Ø19	10997
21	170	21	20	1,7	Ø21	11006
23	170	23	22	1,7	Ø23	11013

The size of the tube is specified by the valve Ø.

Code Catalog 600216  
Rev.01 - 07/2018