


Organic Tubular Arterial Graft Valvulated Tube 300 Corrugated



Experience
Quality
Performance


 **BRAILE**[®]
BIOMÉDICA

Description:

The Organic Tubular Arterial Graft Valvulated Tube 300 Corrugated is made with a bovine pericardium, molded in a special support, treated with glutaraldehyde and conserved in formaldehyde 4%. This treatment gives it proper characteristics of resistance, flexibility and lack of antigenicity.

Indication:

Organic Tubular Arterial Graft Valvulated Tube 300 Corrugated is indicated for corrections of aneurysms of the ascending aorta with valvar insufficiency.



TECHNICAL SPECIFICATIONS (mm)						
Size	A ± 10	ØB ± 0,5	ØC ± 1,0	D ± 0,5	Valve	Code
25	120	27	25	1,7	Ø25	11556
27	120	29	27	1,7	Ø27	11563
29	120	31	29	1,7	Ø29	11570
31	120	33	31	1,7	Ø31	11587
33	120	35	33	1,7	Ø33	11594

The size of the tube is specified by the valve Ø.

Code Catalog 600222
Rev.01 - 07/2018

