

Organic Tubular Arterial Graft Valvulated Tube 300 Smooth



Experience
Quality
Performance


 **BRAILE**[®]
BIOMÉDICA

Description:

The Organic Tubular Arterial Graft Valvulated Tube 300 Smooth is made of bovine pericardium treated with glutaraldehyde and conserved in formaldehyde 4%. This treatment gives it proper characteristics of resistance, flexibility and lack of antigenicity.

Indication:

Organic Tubular Arterial Graft Valvulated Tube 300 Smooth is indicated for corrections of aneurysms of the ascending aorta with valvar insufficiency



TECHNICAL SPECIFICATIONS (mm)					
Size	A ± 10	ØB ± 0,5	ØC ± 1,0	Valve	Code
25	110	27	25	Ø25	10935
27	110	29	27	Ø27	10942
29	110	31	29	Ø29	10959
31	110	33	31	Ø31	10966
33	110	35	33	Ø33	10973

The size of the tube is specified by the valve Ø.

Catalogue Code 600220
Rev. 01 - 07/2018

