

# Organic Tubular Arterial Graft Simple Smooth Tube



Experience  
Quality  
Performance


 **BRAILE**<sup>®</sup>  
BIOMÉDICA

**Description:**

The Organic Tubular Arterial Graft Simple Smooth Tube is made of bovine pericardium treated with glutaraldehyde and conserved in formaldehyde 4%. This treatment gives it proper characteristics of resistance, flexibility and lack of antigenicity.

**Indication:**

The Organic Tubular Arterial Graft in Bovine Pericardium Simple Smooth Tube is indicated for the reconstruction of aneurysms or aneurysmatic diseases, aortoiliac, aortobifemoral revascularization, among others.

**TECHNICAL SPECIFICATIONS (mm)**

Size	A ± 10	ØB ± 0,5	ØC ± 0,5	Code
19	120	19	18	7476
21	120	21	20	7517
23	120	23	22	7559
25	120	25	24	8051
27	120	27	26	8064
29	120	29	28	8077
31	120	31	30	8080
33	120	33	32	8093