

# Bovine Pericardium Patch Pledges




Experience  
Quality  
Performance

### Description:

The Bovine Pericardium Patch Pledges is made with bovine pericardium previously treated with glutaraldehyde and conserved in formaldehyde 4%. This treatment gives it proper characteristics of resistance, flexibility and lack of antigenicity.

### Indication:

The Bovine Pericardium Patch Pledges is used in making and reinforcing the patch for use in extracorporeal circulation, miocardic lesions with an haemostatic purpose, coronary cerclage in procedures without extracorporeal circulation and surgical corrections of the aorta diseases.



TECHNICAL SPECIFICATIONS (cm)			
Size	A $\pm 1$	B $\pm 1$	Code
1 x 0,5	1,0	0,5	4446

Polyacetal Support with 10 units.