

Kit for Autotransfusion



www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA

Description: The Kit for Autotransfusion comprises a reservoir made of polycarbonate, with a 2,000 ml-capacity. It is provided with a system for bubble rupture and a filter made of polyester for retaining particles with a diameter greater than 150µm and PVC tubes.

Indication: The Kit for Autotransfusion is intended for surgical procedures or post-surgery for collecting extravasated blood. This blood will be filtered for later patient reinfusion. This technique uses vacuum for aspirating bleeding.



Code

225690

225690

Kit for Autotransfusion HPO

225645

Kit for Autotransfusion

214944

Kit for Autotransfusion BPSP

225672

Kit for Autotransfusion SPA

Technical Specifications

Raw Material	Polycarbonate
Maximum Volume of Reservoir	2000ml
Maximum Volume of Reservoir	250ml
Connector para Aspirators	1/4"
Aspirated Blood Filter	150 µm
Reservoir Blood Outlet	1/4"



Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I

Zip Code 15091-450 • São José do Rio Preto • SP • Brazil

Phone #55 17 2136 7000 • Fax #55 17 2136 7040

Federal Tax ID: 52.828.936/0001-09

Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010

ANVISA Reg. nº 10159030028

www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA