


# Gas Mixer MGB-1

Simplicity

Accuracy

Reliability



[www.braile.com.br](http://www.braile.com.br)

 **BRAILE**<sup>®</sup>  
BIOMÉDICA




**Description:** The Gas Mixer MGB-1 is a device which allows the accurate control of the mix of medicinal compressed air and oxygen, allowing concentration adjustments from 21 to 100%.

**Indication:** The Gas Mixer MGB-1 is indicated for use with membrane oxygenators in surgical procedures with cardiopulmonary bypass.

Code

44561

Gas Mixer MGB-1



Male air connector (standard)	3/4" - 16 fpp
Male oxygen connector (standard)	9/16" - 18 fpp
Width	30 cm
Length	15 cm
Height	18 cm
Weight	1,2 kg
Oxygen supply pressure	50 psi (345 kPa)
Compressed air supply pressure	50 psi (345 kPa)
Concentration	21 a 100%
Accuracy	5% for flow ranging from 0 to 81/min

## Accessories

Flow meter Connection 1/4 1 to 7 L/min
Compressed air hose Length 5 m (±0.1m)
Oxygen hose Length 5 m (±0.1m)



Catalogue Code 600323  
Rev.

### Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarrá I  
Zip Code 15091-450 • São José do Rio Preto • SP • Brazil  
Phone #55 17 2136 7000 • Fax #55 17 2136 7040  
Federal Tax ID: 52.828.936/0001-09  
Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010  
ANVISA Reg. nº 10159030046  
www.braile.com.br

**BRAILE**<sup>®</sup>  
BIOMÉDICA