

# Digital Clinical Thermometer AG-2000



[www.braile.com.br](http://www.braile.com.br)

 **BRAILE**<sup>®</sup>  
BIOMÉDICA

The Digital Clinical Thermometer AG-2000 is composed of an electronic circuit, visor and supply battery. It is indicated to measuring physiological temperatures. This thermometer obtains its readings via the proper Braille Biomédica sensors.

## Code

117140  
Digital Clinical  
Thermometer AG-2000

### Technical Specifications

Supply	9 VDC (alkaline battery)
Scale	0 - 70°C
Resolution	0,1°C
Accuracy	± 0,4°C
Dimensions	163x105x90mm including clamp
Weight	374g including clamp



#### Braille Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I

Zip Code 15091-450 • São José do Rio Preto • SP • Brazil

Phone #55 17 2136 7000 • Fax #55 17 2136 7040

Federal Tax ID: 52.828.936/0001-09

Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010

ANVISA Reg. nº 10159030029

[www.braile.com.br](http://www.braile.com.br)

 **BRAILE**<sup>®</sup>  
BIOMÉDICA