

# Temperature Sensor



[www.braile.com.br](http://www.braile.com.br)

 **BRAILE**<sup>®</sup>  
BIOMÉDICA

The temperature sensors are internally shielded thermoplastic cables, with a five tract connector with gold terminals and a small-sized PT-100 type pad, enabling faster measurements.

## Technical Specifications

Scale	0 - 70°C
Connection	Redel 5 Vias
Accuracy	± 0,4°C
Length	3m
Sterilizable in Ethylene Oxide	

Code

117168 Temperature Sensor



Catalog Code 600329  
Rev.

### Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I  
Zip Code 15091-450 • São José do Rio Preto • SP • Brazil

Phone #55 17 2136 7000 • Fax #55 17 2136 7040

Federal Tax ID: 52.828.936/0001-09

Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010

ANVISA Reg. nº 10159030075

[www.braile.com.br](http://www.braile.com.br)

 **BRAILE**<sup>®</sup>  
BIOMÉDICA