

WIRED ARTERIAL CANNULA



Thin Wall
Transparent and Flexible Body
Accurate Positioning in the Artery
Kink-Resistant

Description: The Wired Arterial Cannula is made of transparent polyvinyl chloride (PVC), comprising a tip, body and distal end. The tip is provided in straight and curved shapes; it may have a flange for aiding the attachment. The body has a flexible wall with a stainless steel wire reinforcement which avoids kinks while handling it, besides the mark for an accurate positioning in the artery. The distal end may have a connector with a luer lock adapter for removal of samples and measurement of blood pressure. Each package contains a sterile and pyrogen-free unit.

Indication: The Wired Arterial Cannula is intended for cannulation of arteries to the bypass cardiopulmonary perfusion system, providing the patient return of blood.

Description	Diameter	Lenght	Connection	Connector	Code	Situation
Wired Arterial Cannula (Straight Tip) With Flange	18 FR	30 cm	3/8"	3/8" x 3/8" com luer	197018	Estandard Product
Wired Arterial Cannula (Straight Tip) With Flange	20 FR	30 cm	3/8"	3/8" x 3/8" com luer	197027	Estandard Product
Wired Arterial Cannula (Straight Tip) With Flange	22 FR	30 cm	3/8"	3/8" x 3/8" com luer	197045	Estandard Product
Wired Arterial Cannula (Straight Tip) With Flange	24 FR	30 cm	3/8"	3/8" x 3/8" com luer	197054	Estandard Product
Wired Arterial Cannula (Curved Tip) 45°	18 FR	30 cm	3/8"	3/8" x 3/8" com luer	196972	Estandard Product
Wired Arterial Cannula (Curved Tip) 45°	20 FR	30 cm	3/8"	3/8" x 3/8" com luer	196981	Estandard Product
Wired Arterial Cannula (Curved Tip) 45°	22 FR	30 cm	3/8"	3/8" x 3/8" com luer	196990	Estandard Product
Wired Arterial Cannula (Curved Tip) 45°	24 FR	30 cm	3/8"	3/8" x 3/8" com luer	197009	Estandard Product
Wired Arterial Cannula (Curved Tip) 80°	18 FR	30 cm	3/8"	3/8" x 3/8" com luer	192509	Estandard Product
Wired Arterial Cannula (Curved Tip) 80°	20 FR	30 cm	3/8"	3/8" x 3/8" com luer	192387	Estandard Product
Wired Arterial Cannula (Curved Tip) 80°	22 FR	30 cm	3/8"	3/8" x 3/8" com luer	192396	Estandard Product
Wired Arterial Cannula (Curved Tip) 80°	24 FR	30 cm	3/8"	3/8" x 3/8" com luer	192402	Estandard Product
Wired Arterial Cannula (Straight Short Tip) With Flange	08 FR	13 cm	1/4"	—	155876	Made to order
Wired Arterial Cannula (Straight Short Tip) With Flange	10 FR	13 cm	1/4"	—	155885	Made to order
Wired Arterial Cannula (Straight Short Tip) With Flange	12 FR	13 cm	1/4"	—	155894	Made to order
Wired Arterial Cannula (Straight Short Tip) With Flange	14 FR	15 cm	1/4"	—	155900	Made to order
Wired Arterial Cannula (Straight Short Tip) With Flange	16 FR	15 cm	1/4"	—	155919	Made to order
Wired Arterial Cannula (Straight Short Tip) With Flange	18 FR	15 cm	1/4"	—	155928	Made to order
Wired Arterial Cannula (Straight Short Tip) With Flange	20 FR	15 cm	1/4"	—	155937	Made to order
Wired Arterial Cannula (Straight Short Tip) Without Flange	08 FR	13 cm	1/4"	—	157225	Made to order
Wired Arterial Cannula (Straight Short Tip) Without Flange	10 FR	13 cm	1/4"	—	157234	Made to order
Wired Arterial Cannula (Straight Short Tip) Without Flange	12 FR	13 cm	1/4"	—	157243	Made to order
Wired Arterial Cannula (Straight Short Tip) Without Flange	14 FR	15 cm	1/4"	—	157261	Made to order
Wired Arterial Cannula (Straight Short Tip) Without Flange	16 FR	15 cm	1/4"	—	157270	Made to order
Wired Arterial Cannula (Straight Short Tip) Without Flange	18 FR	15 cm	1/4"	—	157289	Made to order
Wired Arterial Cannula (Straight Short Tip) Without Flange	20 FR	15 cm	1/4"	—	157298	Made to order

Description	Diameter	Lenght	Con.	Lenght	Conector	Code	Situation
Wired Arterial Cannula (Elongated Straight Tip)	06 FR	23cm	1/4"	28cm	1/4" x 1/4" com luer	238874	Made to order
Wired Arterial Cannula (Elongated Straight Tip)	08 FR	25cm	1/4"	30cm	1/4" x 1/4" com luer	232283	Estandard Product
Wired Arterial Cannula (Elongated Straight Tip)	10 FR	25cm	1/4"	30cm	1/4" x 1/4" com luer	232292	Estandard Product
Wired Arterial Cannula (Elongated Straight Tip)	12 FR	25cm	1/4"	30cm	1/4" x 1/4" com luer	232308	Estandard Product
Wired Arterial Cannula (Elongated Straight Tip)	14 FR	25cm	1/4"	30cm	1/4" x 1/4" com luer	232317	Estandard Product
Wired Arterial Cannula (Elongated Straight Tip)	16 FR	25cm	1/4"	30cm	1/4" x 1/4" com luer	232326	Estandard Product

Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarrá I
 Zipe Code 15091-450 • São José do Rio Preto • SP • Brazil
 Phone #55 17 2136 7000 • Fax #55 17 2136 7040
 Federal Tax ID: 52.828.936/0001-09
 Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010
 ANVISA Reg. nº 10159030070
 www.braille.com.br

