

NON-WIRED ARTERIAL CANNULA



Thin Wall
Transparent and Flexible Body
Accurate Positioning in the Artery

Description: The Non-wired Arterial Cannula is made of transparent polyvinyl chloride (PVC), comprising a tip, body and distal end. The tip is provided in straight and curved shapes; it may have a flange for aiding the attachment. The body has a thin and flexible wall with marks for an accurate positioning in the artery. The distal end may have a connector with a luer lock adapter for removal of samples and measurement of blood pressure. Each package contains a sterile and pyrogen-free unit.

Indication: The Non-Wired Arterial Cannula is intended for cannulation of arteries to the bypass cardiopulmonary perfusion system, providing the return of blood to the patient.



Catalogue Code 600366
Rev. 00

Description	Diameter	Length	Connection	Connector	Code	Situation
Non-wired Arterial Cannula (Straight Tip)	08 FR	16 FR	1/4"	-	235778	Estandar Product
Non-wired Arterial Cannula (Straight Tip)	10 FR	16 FR	1/4"	-	235787	Made to Order
Non-wired Arterial Cannula (Straight Tip)	12 FR	16 FR	1/4"	-	235796	Made to Order
Non-wired Arterial Cannula (Straight Tip)	14 FR	16 FR	1/4"	-	235802	Made to Order
Non-wired Arterial Cannula (Straight Tip)	16 FR	16 FR	1/4"	-	235811	Made to Order
Non-wired Arterial Cannula (Straight Tip)	16 FR	16 FR	3/8"	-	235820	Made to Order
Non-wired Arterial Cannula (Straight Tip)	18 FR	16 FR	3/8"	-	235839	Estandar Product
Non-wired Arterial Cannula (Straight Tip)	20 FR	16 FR	3/8"	-	235848	Made to Order
Non-wired Arterial Cannula (Straight Tip)	22 FR	16 FR	3/8"	-	235857	Estandar Product
Non-wired Arterial Cannula (Straight Tip)	24 FR	16 FR	3/8"	-	235866	Made to Order
Non-wired Arterial Cannula (Curved Tip) Without Flange	20 FR	18 FR	3/8"	3/8" x 3/8" com luer	210720	Made to Order
Non-wired Arterial Cannula (Curved Tip) Without Flange	22 FR	18 FR	3/8"	3/8" x 3/8" com luer	210739	Made to Order
Non-wired Arterial Cannula (Curved Tip) Without Flange	24 FR	18 FR	3/8"	3/8" x 3/8" com luer	210748	Made to Order
Non-wired Arterial Cannula (Curved Tip) With Flange	20 FR	18 FR	3/8"	3/8" x 3/8" com luer	132851	Made to Order
Non-wired Arterial Cannula (Curved Tip) With Flange	22 FR	18 FR	3/8"	3/8" x 3/8" com luer	132842	Estandar Product
Non-wired Arterial Cannula (Curved Tip) With Flange	24 FR	18 FR	3/8"	3/8" x 3/8" com luer	131432	Made to Order

Braille Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
 Zip Code 15091-450 • São José do Rio Preto • SP • Brazil
 Phone #55 17 2136 7000 • Fax #55 17 2136 7040
 Federal Tax ID: 52.828.936/0001-09
 Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010
 ANVISA Reg. nº 10159030071
www.braille.com.br

