

Non-wired Venous Cannula

Transparent
Flexible



www.braile.com.br

 **BRAILE**[®]
BIOMÉDICA

Description: The Non-wired Venous Cannula is made of transparent and flexible PVC, which is suitable for medical and pharmaceutical use and sterilizable in ethylene oxide (ETO). It is provided with an angled tip, beveled shape, with a side hole, cylindrical body and distal end with a proper design for coupling the standard measure connectors.

Indication: The Non-wired Venous Cannula is intended for venous drainage in surgical procedures performed with cardiopulmonary bypass.



Description	Diameter	Length	Connection	Connector	Code
Non-wired Venous Cannula Tip 45° 3/8"	12,8mm	30cm	3/8"	3/8" x 3/8"	301387
Non-wired Venous Cannula Tip 90° 12Fr	3,9mm	25cm	1/4"	-	167561
Non-wired Venous Cannula Tip 90° 13Fr	4,3mm	27cm	1/4"	-	167570
Non-wired Venous Cannula Tip 90° 14Fr	4,6mm	28cm	1/4"	-	167589
Non-wired Venous Cannula Tip 90° 15Fr	5,0mm	23cm	1/4"	-	167598
Non-wired Venous Cannula Tip 90° 17Fr	5,6mm	28cm	1/4"	-	167604
Non-wired Venous Cannula Tip 90° 20Fr	6,7mm	26cm	1/4"	-	167613
Non-wired Venous Cannula Tip 90° 22Fr	7,3mm	32cm	1/4"	-	167622
Non-wired Venous Cannula Tip 90° 30Fr	9,5mm	32cm	1/4"	-	310453

Catalogue Code 601211
Rev.00

Braile Biomédica Indústria, Comércio e Representações Ltda.

Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
 Zipe Code 15091-450 • São José do Rio Preto • SP • Brazil
 Phone #55 17 2136 7000 • Fax #55 17 2136 7040
 Federal Tax ID: 52.828.936/0001-09
 Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010
 ANVISA Reg. nº 10159030090
www.braile.com.br

