


SURGICAL STENT-GRAFT

Self-expandible Surgical Stent-Graft Endoprosthesis




Self-expansible metallic endoprosthesis with low porosity polyester coating and Nitinol structure.
 For procedures with circulatory arrest (CEC).
 Short release catheter for easy handling in the surgical field.




Suture ◀

Sutures in polyester for the fixation of the metallic skeleton




Coating ◀

Low porosity tubular polyester



400mm



B

A	B
Diameter 26, 28, 30, 32 e 34 mm	Length 75 - 115 - 150 mm

***Ask about availability of other dimensions.**