

Digital Triple Chronometer C3000



Simplicity
Accuracy
Reliability

The Digital Triple Chronometer c3000 Braile Biomedica comprises three high accuracy chronometer coupled to a single box and having independent operation. It works with a 9.6 V built-in battery which can be recharged and/or connected to the power grid, providing a continuous operation even when the local network power source is not available; microprocessed circuit for better reliability and clamp for attachment with two positions (vertical and horizontal). The Chronometer is intended for recording the time of the many surgical steps, such as: time of surgery, time of campling, time of cardioplegia, etc.



Code

353856
Digital Triple Chronometer

Technical Specifications

Precisão	99,0% (25°C)
Display capacitu	99:59:59 (h : m :s)
Duration of Battery	2 hours
Alimentation	power grid or battery

Catalogue Code 600317
Rev. 01

Braile Biomédica Indústria, Comércio e Representações Ltda.
Av. Pres. Juscelino K. de Oliveira, 1505 - Jardim Tarraf I
CEP 15091-450 - São José do Rio Preto - SP - Brasil
Telefone 17 2136 7000 - Fax 17 2136 7040
Other countries - Phone # (55) 17 2136 7000 - Fax 17 2136 7040
CNPJ: 52.828.936-0001-09
Responsável Técnico: Vladimir D. A. Ramirez - CRF/SP 9010
www.braile.com.br

