

# SUCTION DRAIN



**Description:** The Suction Drain comprises a tip coupled to a flexible extension, and a vacuum adapter in the opposite end. All components are made of materials which are suitable for medical and pharmaceutical use. The product works connected to a vacuum source, resulting in a suction effect that permits aspiration and drainage of fluids by means of the tip.

**Indication:** The Suction Drain is intended for draining fluids (blood, secretions or solutions), by aspiration during surgical procedures.

## Technical Specifications

	<b>AD 300</b>
Tube Material	PVC
Tube Length	300 cm
Tip Material	PVC
Vacuum Adapter Material	Silicon

\* See image

## Code

**11037**

Suction Drain - AD 300

