

# RECIRCULATION FILTER




**Description:** The Recirculation Filter is composed by a plastic frame of polymethylmethacrylate (PMMA), with a filtering membrane of polyvinylidene fluoride (PVDF) with a 5µm porosity in its interior. All components are appropriate for medical-pharmaceutical applications, and they are sterilizable in Ethylene Oxide.

**Indication:** The Recirculation Filter is indicated for perfusate filtering during the priming procedure of the extracorporeal blood circulation circuit, in cardiovascular surgeries.

<b>Raw Material</b>	Polymethylmethacrylate (PMMA)
<b>Filtering Membrane</b>	Polyvinylidene Fluoride
<b>Connection</b>	3/8" ou 1/2"
<b>Membrane Porosity</b>	5 µm

**Code**

3478  
Recirculation Filter



Código Catálogo 600269  
Rev. 01

Braile Biomédica Indústria, Comércio e Representações Ltda.  
Av. Pres. Juscelino K. de Oliveira, 1505 - Jardim Tarraf I  
Zip Code 15091-450 - São José do Rio Preto - SP - Brasil  
Phone # (55) 17 2136 7000 - Fax # (55) 17 2136 7040  
Federal Tax ID: 52.828.936-0001-09  
Technician in charge: Vladimir D. A. Ramirez - CRF/SP 9010  
Registro Anvisa nº 10159030087  
www.braile.com.br

