

# KIT OF CANNULAS CEC



**Description:** The Kit of Cannulas CEC comprises a set of cannulas manufactured in extruded, transparent and flexible polyvinyl chloride (PVC), placed on a plastic case (support).

**Indication:** The Kit of Cannulas CEC is indicated to promote the cannulation in surgical procedures conducted with cardiopulmonary bypass (CBP).

## Technical Specifications

Model: ADULT I - Code 3829			
Cannulas	Diameter [inches]	Length [cm]	Quantity [units]
Drainage 3/8"	3/8"	30,0 cm	2
Drainage 3/8" (two section)	3/8"	32,0 cm	1
Suction 1/4"	1/4"	37,0 cm	1
Accessories			
Conducting Cannula	3/8"	19,0 cm	1
Conducting Cannula	1/4"	19,0 cm	1
Straight Connector	1/4"x1/4"	-	1
Flat Lace	-	140 cm	1
Clamp Protector (De Bakey)	-	2 cm	4
Clamp Protector ( Kelly/Halsted (mosquito)	-	1 cm	10
Humel Plug	-	-	5
Plastic Thread Guide	-	-	5
Cord	-	120 cm	3
Metallic Thread Guide	-	-	3
Male Double Luer Connector	-	-	3

Model: ADULT II - Code 242352			
Cannulas	Diameter [inches]	Length [cm]	Quantity [units]
Drainage 3/8"	3/8"	30,0 cm	2
Drainage 3/8" (duplo estágio)	3/8"	32,0 cm	1
Suction 1/4"	1/4"	37,0 cm	1
Accessories			
Straight Connector	1/4"x1/4"	-	1
Needle	-	-	1
Jelco	14 G	-	1
Flat Lace	-	80 cm	4
Clamp Protector (De Bakey)	-	2 cm	10
Clamp Protector ( Kelly/Halsted (mosquito)	-	1 cm	10
Humel Plug	-	-	5
Plastic Thread Guide	-	-	5
Cord	-	80 cm	6
Metallic Thread Guide	-	-	3
Male Double Luer Connector	-	-	4
Reperfusor	-	-	1
Connector Y	1/2"x3/8"x3/8"	-	1
Preparation Balls (pipocas)	-	-	10

Model: PEDIATRIC - Code 136826			
Cannulas	Diameter [inches]	Length [cm]	Quantity [units]
Drainage 1/4"	1/4"	34,0 cm	2
Suction 1/4"	1/4"	37,0 cm	1
Accessories			
Straight Connector	1/4"x1/4"	-	1
Conducting Cannula	1/4"	19,0 cm	1
Multi-Hole Drain	1/4"	34,0 cm	1
Flat Lace	-	130 cm	1
Clamp Protector (De Bakey)	-	2 cm	4
Clamp Protector ( Kelly/Halsted (mosquito)	-	1 cm	10
Humel Plug	-	-	5
Plastic Thread Guide	-	-	5
Cord	-	120 cm	3
Male Double Luer Connector	-	-	3

Model: INFANTILE - Code 3803			
Cannulas	Diameter [inches]	Length [cm]	Quantity [units]
Drainage 3/16"	3/16"	34,0 cm	2
Suction 1/4"	1/4"	37,0 cm	1
Accessories			
Straight Connector	1/4"x1/4"	-	1
Conducting Cannula	1/4"	19,0 cm	1
Flat Lace	-	80 cm	1
Clamp Protector (De Bakey)	-	2 cm	4
Clamp Protector ( Kelly/Halsted (mosquito)	-	1 cm	10
Humel Plug	-	-	5
Plastic Thread Guide	-	-	5
Cord	-	120 cm	3
Multi-Hole Drain	1/4"	34,0 cm	1
Male Double Luer Connector	-	-	3

Model: NEONATO - Codes: 10Fr = 136817   12Fr = 184180   14Fr = 184199			
Cânulas	Diâmetro [polegadas]	Comprimento [cm]	Quantidade [unidades]
Drainage 3/16"	3/16"	23,5 cm	3
Suction 3/16"	3/16"	27,0 cm	1
Accessories			
Arterial cannula	10, 12 ou 14 Fr	-	1
Multi-Hole Drain	1/8"	30,0 cm	1
Straight Connector	3/16"x1/4"	-	2
Clamp Protector ( Kelly/Halsted (mosquito)	-	1 cm	10
Humel Plug	-	-	5
Metallic Thread Guide	-	-	1
Silicon Tube	-	10 cm	5
Jelco	14 G	-	1
Jelco	16 G	-	1
Cord	-	50 cm	5