



OxyPrime ECMO with BRCoating®



Código Catálogo nº 618027
Rev. 01 - 21 jul 2021

Braile Biomédica Indústria, Comércio e Representações Ltda.
 Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
 CEP 15091-450 • São José do Rio Preto • SP • Brasil
 Telefone 17 2136 7000 • Fax 17 2136 7040
 Other countries • Phone # (55) 17 2136 7000 • Fax 17 2136 7040
 CNPJ: 52.828.936-0001-09
 Responsável Técnico: Vladimir D. A. Ramirez • CRF/SP 9010
 Registro Anvisa nº 1015903011
 www.braile.com.br



The OxyPrime ECMO is indicated to long term procedures as pulmonar and/or cardiac support.

The polymethylpentene fibers perform an excellent supply of oxygen and an efficient removal of Co2.

It has an optimized heat exchanger with a surface area of 0.3 m2.

BRCoating®

- decrease the interation between blood and circuit
- free of animal origin substance
- make with synthetic albumin



Technical Specifications

House Material	Polycarbonate
Blood flow rate	1,0 – 7,0L/min
Prime Volume	260ml


Gás Exchange

Fiber Material	Polymethylpentene
Efetive membrane area	1,86 m ²
Sensor de temperatura	YSI
Blood outlet connector	3/8"

Heat Exchange

Fiber material	Polyester
Effective membrane area	0,3 m ²
Blood inlet connector	3/8"
Water inlet/outlet connector	Hansen

Oxygen Transfer



CO₂ Removal (1:1)

