



AORTIC

Size	Measures (mm)					Code
	ØA (±0,5)	ØB (±0,5)	ØC (±1,0)	D (±0,5)	E (±0,5)	
17	12	17	21	05	09	610650
19	14	19	23	06	10	610651
21	16	21	25	07	11	610652
23	18	23	27	08	12	610653
25	20	25	29	09	13	610654
27	22	27	31	10	14	610656
29	24	29	33	11	15	610658



MITRAL

Size	Measures (mm)					Code
	ØA (±0,5)	ØB (±0,5)	ØC (±1,0)	D (±0,5)	E (±0,5)	
25	20	25	29	10	13	610655
27	22	27	31	11	14	610657
29	24	29	33	12	15	610659
31	26	31	35	13	16	610660
33	28	33	37	13	16	610661
35	30	35	39	14	17	610662



VIVERE®

BOVINE PERICARDIUM ORGANIC
VALVULAR BIOPROSTHESIS

EXCLUSIVE
REALOG®
ANTI-CALCIFICATION
TREATMENT

PROLONGED
DURABILITY

LOW PROFILE


NEW AORTIC DESIGN
FOR BETTER ADJUSTMENT
OF POSITIONING

To realize a purchase order or to know more about our products, accessories and therapeutic solutions, please contact your sales representative.

- ▶ Polyacetal ring coated with bovine pericardium, providing:
 - ▶ Interface between the ring and leaflets less susceptible to wear, reducing the probability of structural failure
 - ▶ Valvular surface with greater hemocompatibility
- ▶ Leaflets made from a single piece of bovine pericardium, resulting in:
 - ▶ Better wear resistance and functional fatigue
 - ▶ More homogeneous behavior in relation to performance
 - ▶ Guarantee of morphological homogeneity in relation to the distribution of collagen, elastic fibers and mechanical properties
- ▶ Formaldehyde storage
- ▶ **Anti-calcification treatment** by the exclusive **REALOG® process**




AORTIC




REALOG® Process

The process adds a treatment with glutamic acid (Glu), followed by oxidation, resulting in the blocking of free aldehyde groups and in the reduction of cross-links of polymer glutaraldehyde, which reduces the post-implant calcification index, ensuring a longer valve life.



MITRAL



Azeredo LG et al. Late outcome analysis of the Braille Biomedica pericardial valve in the aortic position. Braz J Cardiovasc Surg 2014;29(3):316-21.
 Braille-Sternieri MCVB. Pericárdio bovino tratado com ácido glutâmico para confecção de próteses valvulares cardíacas: estudos *in vitro* e *in vivo*. Thesis, 2013.
 Braille MCVB et al. In vitro properties and performance of glutaraldehyde-crosslinked bovine pericardial bioprostheses treated with glutamic acid. Artif Organs 2011;35(5):497-501.
 Braille DM et al. Aortic valve replacement using bovine pericardial bioprostheses: 12 years of experience. Braz J Cardiovasc Surg 2003;18(3):217-20.
 Braille DM et al. IMC bovine pericardial valve: 11 years. J Cardiac Surg 1991;6(4):580-8.