

# KIT FOR AUTOTRANSFUSION





**Description:** The Autotransfusion Set consists of a polycarbonate reservoir with a capacity of 2,000ml. It has a system for breaking bubbles and a polyester filter to retain particles with a diameter greater than 150µm and PVC tubes.

**Indication:** The product is indicated for surgical or post-operative procedures to collect extravasated blood. This blood will be filtered for subsequent reinfusion to the patient. This technique uses a vacuum to aspirate bleeding.



**Code**

225645  
Kit for Autotransfusion

SPA Autotransfusion System	Autotransfusion System	BPSP Autotransfusion System
Raw material: Polycarbonate	Raw material: Polycarbonate	Raw material: Polycarbonate
Reservoir Capacity: (2000 ml)	Reservoir Capacity: (2000 ml)	Bottle Capacity: (1000 ml)
PVC vacuum tube ¼ (Outlet connector with female luer)	Connector for vacuum cleaners: ¼	Catheter: 3/8 x 45cm
Aspirated Blood Filter: 150um	Aspirated Blood Filter: 150um	Suction drain: AD 300
Exit from the Reservoir: ¼	Exit from the Reservoir: ¼	
Acrylic Vacuum	Connector: Y 3/8 x 3/8 x 3/8	
Connector: Y 3/8 x 3/8 x 3/8		

Código Catálogo 6002E7  
Rev. 02 - Jan/2022