

OxyPrime[®]

Membrane Oxygenator




Low Priming
Volume

New design

More Compact

High Performance



Technique Especification





Material	Polycarbonate
Maximum Blood Flow	7,0l/min
Total Volume Priming	250 ml

Gas Exchanger

Fiber Material	Polypropylene
Effective Membrane Area	1,95 m ²
Temperature Sensor	YSI
Blood Outlet Connector	3/8"
Gas Inlet Connector	1/4"
Recirculation Connector	luer lock
Arterial Sample Connector	luer lock
Air Removal Connector	luer lock

Heat Exchanger

Fiber Material	Polyester
Effective Membrane Area	0,4 m ²
Maximum Water Temperature	42 C
Maximum Water Pressure	40 PSI
Blood Inlet Connector	3/8"
Water Inlet/Outlet Connector	Hansen (1/2")



Commercial Information

OxyPrime 7000

Code 617484

For orders or more information about any of our products, accessories or therapeutic solutions, contact your local sales representative.

Braille Biomédica Indústria, Comércio e Representações Ltda.
 Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
 Zip Code 15091-450 • São José do Rio Preto • SP • Brasil
 Phone # (55) 17 2136 7000 • Fax # (55) 17 2136 7040
 Federal Tax ID: 52.828.936-0001-09
 Technician in charge: Vladimir D. A. Ramirez • CRF/SP 9010
 Registro Anvisa nº 10159030106
 www.braille.com.br

