

KIT OF TUBES FOR ECMO



Description: The Tubing Sets for ECMO with BRCoating is indicated for long-term procedures as pulmonary and/or cardiac support. The tubes are pre assembled tray, made with polyvinyl chloride (PVC), connectors and medical level, non-toxic and biocompatible materials and has the coating treatment BRCoating.



BRCoating

- decrease the interaction between blood and circuit
- free of animal origin substance
- make with synthetic albumin

	Code	Model	Venous Line	Arterial Line	Centrifugal Line
Tubing Set for Adult ECMO with BRCoating	617327	Adult	3/8"	3/8"	3/8"