


AORTIC AND MITRAL VALVE METERS



The heart valve meters are used to determine the diameter of the annulus, defining the ideal size of the prosthesis to be implanted, in order to meet different anatomical characteristics of the patients.

The set is made of 304 stainless steel and contains 06 valve meters and 01 holder cable, the meters are 03 units for the aortic position and 03 units for the mitral position



Sizes (mm)						
Aortic	19	21	23	25	27	29
Mitral	25	27	29	31	33	35

Código 625920
Rev. 00 - Set. 2023