



# OxyPrime ECMO with BRCoating®



Código Catálogo nº 618027 Rev. 01

Braile Biomédica Indústria, Comércio e Representações Ltda.  
 Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I  
 CEP 15091-450 • São José do Rio Preto • SP • Brasil  
 Telefone 17 2136 7000 • Fax 17 2136 7040  
 Other countries • Phone # (55) 17 2136 7000 • Fax 17 2136 7040  
 CNPJ: 52.828.936-0001-09  
 Responsável Técnico: Vladimir D. A. Ramirez • CRF/SP 9010  
 Registro Anvisa nº 1015903011  
 www.braile.com.br



The OxyPrime ECMO is indicated to long term procedures as pulmonar and/or cardiac support.

The polymethylpentene fibers perform an excellent supply of oxygen and an efficient removal of Co2.

It has an optimized heat exchanger with a surface area of 0.3 m2.

## BRCoating®

- decrease the interation between blood and circuit
- free of animal origin substance
- make with synthetic albumin



### Technical Specifications

House Material	Polycarbonate
Blood flow rate	1,0 – 7,0L/min
Prime Volume	260ml


#### Gás Exchange

Fiber Material	Polymethylpentene
Efetive membrane area	1,86 m <sup>2</sup>
Sensor de temperatura	YSI
Blood outlet connector	3/8"

#### Heat Exchange

Fiber material	Polyester
Effective membrane area	0,3 m <sup>2</sup>
Blood inlet connector	3/8"
Water inlet/outlet connector	Hansen

### Oxygen Transfer



### CO<sub>2</sub> Removal (1:1)

