

TEMPERATURE SENSOR



The temperature sensors are internally shielded thermoplastic cables, with a five tract connector with gold terminals and a small-sized PT-100 type pad, enabling faster measurements.

Technical Specifications

| | |
|--------------------------------|--------------|
| Scale | 0 - 70°C |
| Connection | Redel 5 Vias |
| Accuracy | ±0,4°C |
| Length | 3m |
| Sterilizable in Ethylene Oxide | |

Code

117168

Temperature Sensor



Código Catálogo 600329
Rev. 01 - mar 2021

Braille Biomédica Indústria, Comércio e Representações Ltda.
Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I
CEP 15091-450 • São José do Rio Preto • SP • Brasil
Telefone 17 2136 7000 • Fax 17 2136 7040
Other countries • Phone # (55) 17 2136 7000 • Fax 17 2136 7040
CNPJ: 52.828.936-0001-09
Responsável Técnico: Vladimir D. A. Ramirez • CRF/SP 9010
Registro Anvisa nº 10159030075
www.braille.com.br

 **BRAILE**[®]
BIOMÉDICA