

# DIGITAL CLINICAL THERMOMETER AG-2000



The Digital Clinical Thermometer AG-2000 is composed of an electronic circuit, visor and supply battery. It is indicated to measuring physiological temperatures. This thermometer obtains its readings via the proper Braille Biomédica sensors.

Code

117140  
Digital Clinical  
Thermometer AG-2000

#### Technical Specifications

Supply	9 VDC (alkaline battery)
Scale	0 - 70°C
Resolution	0,1°C
Accuracy	±0,4C
Dimensions	163x105x90mm including clamp
Weight	374g including clamp



Código Catálogo 600331  
Rev. 01 - mar 2021

Braille Biomédica Indústria, Comércio e Representações Ltda.  
Av. Pres. Juscelino K. de Oliveira, 1505 • Jardim Tarraf I  
CEP 15091-450 • São José do Rio Preto • SP • Brasil  
Telefone 17 2136 7000 • Fax 17 2136 7040  
Other countries • Phone # (55) 17 2136 7000 • Fax 17 2136 7040  
CNPJ: 52.828.936-0001-09  
Responsável Técnico: Vladimir D. A. Ramirez • CRF/SP 9010  
Registro Anvisa nº 10159030029  
www.braille.com.br

 **BRAILE**<sup>®</sup>  
BIOMÉDICA