

# Vesto<sup>®</sup>




 **BRAILE**

LIFE-SAVING TECHNOLOGIES


The Vesto® Vascular Stent – Endoprosthesis is a device intended for permanent implantation in the peripheral arterial system. It consists of a self-expanding nitinol metallic structure covered with polytetrafluoroethylene (PTFE) and polyester fibers. This structure prevents stent elongation.

Vesto® is indicated for the endovascular treatment of conditions such as occlusions, stenoses, dissections, and aneurysms located in the iliac and superficial femoral arteries.

## Flexible and self-expanding stent



High radiopacity




Gold radiopaque markers  
at the extremities

## Delivery System

**Radiopaque tip** that enables guidance during delivery system navigation and endovascular prosthesis deployment.

**Atraumatic design** that prevents vascular damage during delivery system navigation.



Excellent trackability due to a reinforced, flexible, anti-kink sheath with hydrophilic coating.

## Vesto Long

120 cm delivery system length

Code	Description - Diameter (mm) x Length (mm)	Profile
616898	Stent Vascular - Endoprótese Vesto Longo 6/60	8Fr
616899	Stent Vascular - Endoprótese Vesto Longo 6/80	8Fr
616900	Stent Vascular - Endoprótese Vesto Longo 6/100	8Fr
616901	Stent Vascular - Endoprótese Vesto Longo 6/120	8Fr
616902	Stent Vascular - Endoprótese Vesto Longo 6/140	8Fr
616903	Stent Vascular - Endoprótese Vesto Longo 8/60	9Fr
616904	Stent Vascular - Endoprótese Vesto Longo 8/80	9Fr
616905	Stent Vascular - Endoprótese Vesto Longo 8/100	9Fr
616906	Stent Vascular - Endoprótese Vesto Longo 8/120	9Fr
616907	Stent Vascular - Endoprótese Vesto Longo 8/140	9Fr
616908	Stent Vascular - Endoprótese Vesto Longo 10/60	10Fr
616909	Stent Vascular - Endoprótese Vesto Longo 10/80	10Fr
616910	Stent Vascular - Endoprótese Vesto Longo 10/100	10Fr
616911	Stent Vascular - Endoprótese Vesto Longo 10/120	10Fr
616912	Stent Vascular - Endoprótese Vesto Longo 10/140	10Fr
616913	Stent Vascular - Endoprótese Vesto Longo 12/60	11Fr
616914	Stent Vascular - Endoprótese Vesto Longo 12/80	11Fr
616915	Stent Vascular - Endoprótese Vesto Longo 12/100	11Fr
616916	Stent Vascular - Endoprótese Vesto Longo 12/120	11Fr
616917	Stent Vascular - Endoprótese Vesto Longo 12/140	11Fr
616918	Stent Vascular - Endoprótese Vesto Longo 14/60	12Fr
616919	Stent Vascular - Endoprótese Vesto Longo 14/80	12Fr
616920	Stent Vascular - Endoprótese Vesto Longo 14/100	12Fr
616921	Stent Vascular - Endoprótese Vesto Longo 14/120	12Fr
616922	Stent Vascular - Endoprótese Vesto Longo 14/140	12Fr



## Vesto Short

70 cm delivery system length



Code	Description - Diameter (mm) x Length (mm)	Profile
619267	Stent Vascular - Endoprótese Vesto Curto 6/60	8Fr
619268	Stent Vascular - Endoprótese Vesto Curto 6/80	8Fr
619269	Stent Vascular - Endoprótese Vesto Curto 6/100	8Fr
619270	Stent Vascular - Endoprótese Vesto Curto 6/120	8Fr
619271	Stent Vascular - Endoprótese Vesto Curto 6/140	8Fr
619272	Stent Vascular - Endoprótese Vesto Curto 8/60	9Fr
619273	Stent Vascular - Endoprótese Vesto Curto 8/80	9Fr
619274	Stent Vascular - Endoprótese Vesto Curto 8/100	9Fr
619275	Stent Vascular - Endoprótese Vesto Curto 8/120	9Fr
619276	Stent Vascular - Endoprótese Vesto Curto 8/140	9Fr
619282	Stent Vascular - Endoprótese Vesto Curto 10/60	10Fr
619283	Stent Vascular - Endoprótese Vesto Curto 10/80	10Fr
619284	Stent Vascular - Endoprótese Vesto Curto 10/100	10Fr
619285	Stent Vascular - Endoprótese Vesto Curto 10/120	10Fr
619286	Stent Vascular - Endoprótese Vesto Curto 10/140	10Fr
619297	Stent Vascular - Endoprótese Vesto Curto 12/60	11Fr
619298	Stent Vascular - Endoprótese Vesto Curto 12/80	11Fr
619299	Stent Vascular - Endoprótese Vesto Curto 12/100	11Fr
619300	Stent Vascular - Endoprótese Vesto Curto 12/120	11Fr
619301	Stent Vascular - Endoprótese Vesto Curto 12/140	11Fr
619302	Stent Vascular - Endoprótese Vesto Curto 14/60	12Fr
619303	Stent Vascular - Endoprótese Vesto Curto 14/80	12Fr
619304	Stent Vascular - Endoprótese Vesto Curto 14/100	12Fr
619305	Stent Vascular - Endoprótese Vesto Curto 14/120	12Fr
619306	Stent Vascular - Endoprótese Vesto Curto 14/ 140	12Fr

Follow us on social media:

 /brailebiomedica

 @brailebiomedica



**Braile Biomédica Indústria, Comércio e Representações Ltda.**

Av. Pres. Juscelino K. de Oliveira, 1505 | Jardim Tarraf I

CEP 15091-450 | São José do Rio Preto - SP | Brasil

CNPJ: 52.828.936/0001-09 | Anvisa no. 10159030113

Technical Responsible: Vlademir D. A. Ramirez - CRF-SP 09010